#### USO

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

#### BIBLIOGRAFIA

1. Gotto, AM, Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice, 23; Suppl. 1, 4 (1988).
2. Crouse, JR et al., Studies of low density lipoprotein molecular weight in human with coronary artery disease, J. Lipid Res., 26:566 (1985).
3. Badimon, JJ, Badimon, L., Fuster V., Regression of Atherosclerotic Lesions by HDL. Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 1990; 85:1234-41.
4. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55;767 (1977).
5. Barr, DP, Russ EM, Eder, HA, Protein-lipid relationships in human plasma, Am. J. Med., 11; 480 (1951).
6. Gordon, T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med., 62:707 (1977).
7. Williams, P., et al, HDL and coronary risk factor, Lancet, 1; 72, (1979).
8. Kannel, WB, Castelli, WP, Gordon, T., Cholesterol in the prediction of atherosclerotic disease.; Ann. Intern. Med., 90:85, (1979).
9. National Institutes of Health publication No. 93-3095, Sept. (1993).
10. Special Communication, Executive Summary of 3rd Report of the NCEP. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol (Adult Treatment Panel III), JAMA, Vol. 285, No. 19, May 16, 2001, pages 2486 - 2497.
11. Warnick, GR, Wood, PD, National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary, Clinical Chemistry, Vol. 41, No. 10, 1427-1433 (1995).
12. Kimberly MM, et al. Selection, validation, standardization, and performance of a designated comparison method for HDL-cholesterol for use in the cholesterol reference method laboratory network. Clin Chem 1999; 45:1803-12.
13. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 6, No. 13, August (1986).
14. Young, DS, Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, DC, 1990, 3-104 thru 3-106.
15. Camps, J, Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol, Clinical Chemistry, 1999; 45:685-688.
16. Tietz, NW, Clinical Guide to Laboratory Tests, WB Saunders Co., 1986, p. 256.
17. Carey RN, Garber CC. Evaluation of methods. In: Kaplan LA, Pesce AJ, eds. Clinical Chemistry: theory, analysis and correlation. Third Edition. St. Louis: The CV Mosby Company, St. Louis, MO., 1996: 402-423.
18. Westgard JO, Carey RN, Wold S. Criteria for judging precision and accuracy in method development and evaluation. Clinical Chemistry 1974; 20:825-833.
19. Bachorik PS, Ross JW, for the National Cholesterol Education Program Working Group on Lipoprotein Measurements. National Cholesterol Education Program recommendations for measurement of low-density lipoprotein cholesterol: executive summary, Clin Chem 1995; 41:1414-1433.
20. National Reference system for Cholesterol, CRMLN HDL Cholesterol Protocol, November 2002.
21. Kimberly MM, et al., Selection, validation, standardization, and performance of a designated comparison method for HDL-cholesterol for use in the Cholesterol Reference Method Laboratory Network. Clinical Chemistry 1999; 45:1803-12.

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